

## **Description of the product**

MAMUFACTURER NAME:	DCI INDUSTRY CO.,LTD	
MAMUFACTURER ADDRESS:	26-20 826 beon-gil, Seobu-ro, Jillye-myeon, Gimhae-si, Gyeongsangnam-do, South Korea	
PRODUCT TYPE :	Particle Filtering Half Mask	
PPE REGULATION :	(EU)2016/425-catalog III	
APPICABLE STANDARDS:	EN149:2001+A1:2009	
PERFORMACE CLASSIFICATION:	FFP2 NR ( CE2163)-Certification No : 2163-PPE-1977	
MODEL IDENTIFICATION :	DCI 3D MASK ( DOCTOR SOOMEASY) DCI 3D MASK means that a mask in form of three folds vertically made by DCI INDUSTRY DOCTOR SOOMEASY means that a mask developed by 3 doctors and easy to breathing. (the words " SOOM" means "Breathing" in Korean language ")	
DATE AND REVISION CONTROL :	13 JANUARY 2021 VERSION 0	

#### PRODUCT IMAGE




SUNMITTED BY NAME:	CHOI JONG HWAN	
SUNMITTED BY DATE:	13-Jan-21	


### **Packaging Marking Detail**



### **Design Concept.**



#### Special 4-layer(TSMS) structure mask




Inhaling should be easier than exhaling.


The static electricity of the MB filter should be stored as long as possible.


Since the static electricity is so vulnerable to moisture, the MB filter must be located far away from the mouth.

# **Size Comparison** (Real measured Data)




### **DIFFERENT SIZE MAKE DIFFEERNT VOLUME**





### **Design and Manufacturing Drawing**



## **Test Result of Filtering Efficiency**

### FFP2 grade, but satisfies FFP3 grade

enetration of filter r	naterial: Sodium (	Chloride Testing	No.	
Condition	No. of Sample	Sodium Chloride Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result
(A.R.)	36	0,05		1.
(A.R.)	37	0,14		
(A.R.) 38	0,15	FFP1 ≤ 20 %	Filtering half masks fulfill th	
(S.W.)	1	0,18		requirements of the standar
(S.W.)	2	0,13	FFP2 ≤ 6 %	EN 149:2001 + A1:2009
(S.W.)	3	0,17	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	given in 7.9.2 in range of th FFP1 FFP2 and FFP3
(M.S. T.C.)	10	0,20	FFP3 ≤ 1 %	classes.
(M.S. T.C.)	11	0,34	<del></del>	classes.
(M.S. T.C.)	12	0,29		
onditioning: (M.S.)	Mechanical Streng	th	•	95 L/min = 1,6 dm <sup>3</sup> .sn <sup>-1</sup>
(T.C.)	Temperature Cond	fitioning		SAL CERTIFICA
4 . 70 .	As Received, origi	1-4		CAL

Condition	No. of Sample	Paraffin Oil Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result
(A.R.)	39	0,14		
(A.R.)	40	0,21	]	and a service service and
(A.R.)	A.R.) 41 0,18 EEP1 < 20 %	FFP1 ≤ 20 %	Filtering half masks fulfill the	
(S.W.)	4	0,21		requirements of the standard EN 149:2001 + A1:2009 given in 7.9.2 in range of the FFP1 FFP2 and FFP3 classes.
(S.W.)	5	0,20	FFP2 ≤ 6 %	
(S.W.)	6	0,22	FFP3 ≤ 1 %	
(M.S, T.C.)	13	0,60		
(M.S. T.C.)	14	0,78		
(M.S. T.C.)	15	0,64		
nditioning: (M.S.) Mechanic	al Strength ture Conditioning			

Please refer to the separate test report provided by the manufacturer for other detail test result.

#### **User Information**

#### **DECALARATION**

### Supply of User Information

Declaration Number: DCI-2020108-002

The Manufacture : DCI INDUSTRY CO.,LTD

Company Address: 909-13 Seobu-ro 826beon-gil, Jillye-myeon, Gimhae-si, Gyeongsangnam-do, South Korea. Factory Address: 26-20, Seobu-ro 826beon-gil, Jillye-myeon, Gimhae-si, Gyeongsangnam-do, South Korea.

Provides information about PPE's user information as follows

Product: Particle Filtering Half mask

Product name: DCI 3D MASK ( DOCTOR SOOMEASY MASK)

#### Intended use of the product

This product is intended to be used as a single shift use particle filtering half masks for protection against solid and liquid aerosols.

#### Instructions

- See information supplied by the manufacturer
- Place the particle filtering mask onto the face to cover both nose and chin.
- Hold the particle filtering mask with nose wire section up, cover your nose and mouth.
- Adjust the nose area by pressing the nose wire on the bridge of the nose with your fingers.
- Cover the entire particle filtering mask by both hands and adjust the mask to stick to the face well.

#### Caution

- Single shift use only (NR)
- · Do not use a towel or a facial tissue under the particle filtering mask.
- Do not use if inside of the particle filtering mask is contaminated.
- Do not wash the particle filtering mask.
- Do not deform the shape of the particle filtering mask.
- · Avoid keep the particle filtering mask when wearing it.
- Stop using the particle filtering mask and consult with a doctor in case of any difficulty in breathing while using the particle filtering mask especially for pregnant women, patients with respiratory and cardiovascular disease, children and elderly.
- Do not use in a closed place with less than 18% oxygen concentration.
- As facial hair reduce the protection offered it is recommended that these respirators are used by clean shaven personnel only.
- Do not use for protection against gases.
- Do not remove whilst in contaminated area.

Packing unit: 1pcs

Expiry Date: 36 months after manufacturing date Lot No: Indicated separately

Storage: at room temperature (1~30°C) in dry(maximum relative humidity of storage condition is 70%) ventilated and noncorrosive atmosphere.

Components: Outer Layer, Inner Layer, Filter, Middle Layer(non-woven fabric),

Nose wire (polypropylene coated wire), Ear Loop ( Polyurethane string)

Standard: EN149:2001+A1:2009

Classification: FFP2 NR

C/S: +82-55-902-7110 Email: info@dci-industry.com

This declaration is issued under the responsibility of the manufacture.

2021.01.08

Conforms to EN149:2001+A1:2009 C€2163 FFP2 NR

Product conforms to the requirements of Regulation(EU) 2016/425 of European Parliament of the Council as Personal Protective Equipment (PPE) EU Type Examination(Module B and Module C2) certificate issued by Universal Certification and Surveillance Services Trade Co. in Turkey [N.B No.2163]

DCI INDUSTRY CO., LTD CT INDUSTRY Co., Ltd Jong-Hwan Choi



NB 2163

### **EU TYPE EXAMINATION CERTIFICATE**

Certificate No: 2163-PPE-1977

Respiratory protective devices, filtering half masks to protect against particles manufactured by

#### DCI INDUSTRY CO., LTD.

909-13, Seobu-ro, Jillye-myeon, Gimhae-si, Gyeongsangnam-do, South Korea

are tested and evaluated according to

#### EN 149:2001 + A1:2009 Respiratory Protective Devices -Filtering Half Masks to Protect Against Particles -Requirements, Testing, Marking

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

#### **Product Definition**

Single use particle filtering half mask for protection against solid and liquid aerosols, is a folding type, 4 layered, without valve, ear straps and adjustable nose bar.

Model: DCI 3D MASK Classification: FFP2 NR

For more details, refer technical evaluation report provided to the manufacturer, dated 04.02.2021 and number 2163-KKD-1977.

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.
- Ongoing successful performance in fulfilment of the requirements set out in Personal Protective Equipment Regulation (EU) 2016/425 and harmonised standards, ensured by assessments based on Annex 7 (Module C2) or Annex 8 (Module D) of the regulation.

This certificate is initially issued on 04/02/2021 and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.

CE 2163

Suat KACMAZ
UNIVERSAL CERTIFICATION
Director

validity with the OR



### Splash Test Result. 1/2





Test Report

No.T32120240645SN

Date: Apr 16, 2021

Page 1 of 2

DCI INDUSTRY CO., LTD.

909-13, SEOBU-RO, JILLYE-MYEON, GIMHAE-SI, GYEONGSANGNAM-DO, KOREA

The following samples were submitted and identified by/on behalf of the client as:

DCI 3D MASK (DOCTOR SOOMEASY)

Case No. : CA321202439609 SGS Reference No. : AYUC21-00013 Lot No. / Batch Code : NOT PROVIDED Sample Description : WHITE MASK

: DCI INDUSTRY CO., LTD. Manfacturer

Quantity Submitted : 100 PCS Country of Origin : SOUTH KOREA

Sample Receiving Date : MAR 26, 2021

Testing Period : MAR 26, 2021 - APR 16, 2021

Test Requested	Conclusion
Resistance to penetration by synthetic blood (ASTM F1862/F1862M-17)	PASS

\*\*\*\*\*\*\* FOR FURTHER DETAILS, PLEASE REFER TO THE FOLLOWING PAGE(S) \*\*\*\*\*\*\*\*

Signed for and on behalf of SGS Hong Kong Ltd.

Au Kam Chi, Gigi Technical Manager



### Splash Test Result. 2/2





Test Report No.T32120240645SN Date: Apr 16, 2021 Page 2 of 2

#### Test Results:

Resistance to penetration by synthetic blood ASTM F1862/F1862M-17

Test Side : Outside

Pre-Conditioning : Minimum of 4 hours at 21±5°C and 85±5% R.H.

Test Condition 21±59C and 85±10% R.H.

Test Pressure : 120 mmHg
No of Test Specimen Tested : 32 → LEVEL 2


No of Test Specimen Passed : 31

Test Specimen #	Synthetic Blood Penetration	
1-30, 32	None Seen	
31	Yes	

#### Note:

An acceptable quality limit of 4,0 % is met for a single sampling plan when 29 or more of the 32 tested specimens show "pass" results.

#### Sample Photo:



SGS authenticate the photo on original report only

\*\*\* End of Report \*\*\*

This document is issued by the Company subject to its General Conditions for Service printed overlead, available on request or accessible at <a href="http://www.sgc.com/en/Terms-and-Conditions tor Sectoric Exercises and Conditions For Sectoric Exercises For Sector

Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.

#### **FDA Certification**



#### **FDA Certification**



#### **ISO13485 Certification**

## **N**TREE

Certificate of Registration

### DCI INDUSTRY Co., Ltd.

909-13, Seobu-ro, Jillye-myeon, Gimhae-si, Gyeongsangnam-do, Republic of Korea

has been approved by NTREECERT Co., Ltd. to the following medical devices quality management system standards:

ISO 13485:2016

The scope of medical devices quality management system is applicable to:

Development, Manufacturing and Sales of Dental Mask

#### Valid Date

From 03 November 2020 To 02 November 2023
Former Certificate: 03 November 2020
Certificate No.: NEW 0020 Pt. 1

Certificate No.: NTM-0039 Date of last Issue: 03 November 2020

Rev. 00

nein)

President of NTREE Cert



www.ntreecert.kr NTREE Certification Co.,Ltd



30, Pajangcheon-ro 44beon-gil, Jangan-gu, Suwon-si, Gyeonggi-do, Republic of Korea

### **Certificate of Business Registration**




발 급 번 호 Issuance number 7594-980-6725-864	사 업 자 등 록 증 명 Certificate of Business Regi ( 법인사업자 ) ( Corporate Taxpayer	stration	처리기간 Processing time 즉 시 Immediately
상호 (법 인 명) Name of company	주식회사 디세아이인데스트리 DCI INDUSTRY CO.,LTD	,	
사 업 자 등 목 번 호 Business registration number	297-86-01927		
성 명 (대 표 자) Name of representative	최종환 CHOI JONG HWAN		
주 민 ( 법 인 ) 등 목 번 호 Resident(Corporation) registration number	195511-0246235		
사 업 장 소 계 지 Business Address	경상남도 김해시 전해면 서부로 909-13 908-13, Seobu-ro, Jillye-myeon, Gimhae-si, Gyeongsangnam-do, Republic of Korea		blic of Korea
게 열 일 Date of business commencement	2020년(Year) 07월(Month) 01월(Day)		
사 업 자 등 목 일 Date of business registration	2020년(Year) 07월(Month) 01월(Day)		
업 Business type			
表 Business item	마스크 Face Mask 기계부품 Wholesale of electrical machinery and related materials		
	성명(법인평) Name(Name of company)	주민(A Resident(Busin	업자)등록변호  ess) registration No.
공 등 사 열 자 Joint business owner	채당사항 없습니다 (No Data)	Y	
	위와 같이 증명합니다.		

위와 같이 증명합니다.

I certify that above information is true and correct to the best of my knowledge and belief, 꽃 위 내용은 발급일 현계 상황으로서 추후 변경될 수 있습니다. This information is true as of the issuance date of this certification and but maybe subject to change in the future.


접 수 번 호 Receipt No.	502082209299
당 당 부 서	면원봉사실
Department	Taspayer Service Center
답 당 차	김미정
Staff in Charge	Kim Mi jeong
연 탁 처 Telephone No.	+082-55-320-6224

2020 년 11 월 3 일 Year Month Day



Head of Gimhae District Tax Office









- 본 증명의 위·변조 여부는 발급일로부터 90일 이내「국제청 홈팩스(www.hometax.go.kr) 또는 모바일 홈팩스 > 민원증명(증명발급) > 민원증명 원본확인」에서 발급번호로 확인, 또는 문서 하단의 바코드로 확인이 가능합니다. (공문서를 위·벤조하거나 행사한 자는 10년 이하의 정역에 처할 수 있습니다.)
- \* 본 증명은 홈펙스(www.hometax.go.kr)에서 대면 온라인 서비스를 통해 발급된 증명서입니다.

## Thanks!!

- Our video is provided separately-

