



BLOOD COLLECTION



sol-m.com | info.eu@sol-m.com

EMEA-DIA-382-1022

©2023 Sol-Millennium Europe Sp. z o.o. All rights reserved.









SOL-MILLENNIUM is a vertically integrated global medical supply company group with locations across the globe. We are one of the world's leading manufacturers of needles and syringes and also distribute a comprehensive line of medical devices and supplies.

Our innovative technology reduces needlestick injuries, saves on medication expenses and enhances both clinician and patient experience. Our international sales team is dedicated to delivering knowledge and outstanding customer service. Our goal is to work together as partners in the healthcare industry, focusing on a healthier tomorrow.



- Global Main Office
- Research & Development
- Sales, Marketing & Distribution
- Manufacturing

Table of Contents

	3-6		Safety Blood Collection Tube Holder Set
	7-10		Multi-sample Needle
	11-14		Safety Blood Collection Needle
	15-18		Accessories

This material is a general information resource aimed at demonstrating our offering solely for Healthcare Professionals based in Poland and in some of the EU countries. Its content is directed at and restricted to a professional audience only. It is not a binding offer nor country-specific information, and as the case may be, the products' availability and our offering may vary in the country where you are located. Please contact your Sol Millennium representative if you have questions about the availability of specific products in your area.

A Healthcare Professional must always rely on his or her own professional training, clinical judgment and the patient's condition when deciding whether to use a particular product. This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use.


Unless otherwise indicated, this material is our proprietary property and all copyrights as well as photos, graphics, text, trademarks, service marks, and logos are either owned or licensed by us. This material may not be spread through the Internet, press or social media, as well as reproduced in any form, in whole or in part without our consent.

Sol-Millennium reserves the right to add, change or discontinue products during the validity of this catalogue without previous notice.

This material is distributed by Sol-Millennium Europe Sp z o.o with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIth Economic Department KRS 0000434783, NIP: 5252535720, EU VAT: PL5252535720, Share capital: 500.000 PLN.

BLOOD COLLECTION

Sol-Care™ Safety Blood Collection Tube Holder Set



Sol-Care™ Safety Blood Collection Tube Holder Set is a high-quality and high-tech blood collection device, designed to provide a safe and reliable method for facilitating blood withdrawal from a patient into evacuated blood collection tubes without exposing the healthcare worker to an accidental needle stick injury.

Sol-Care™

● Safety Blood Collection Tube Holder Set

Minimum training requirement

“Using the yellow activation tube is like adding a blood collection tube”.

In-vein retraction

Effectively reduces the risk of needlestick injuries and blood exposure.




Compatible with standard multi-sample needles and blood collection vacuum tubes



Blood Collection Tube Holder



Vacuum Activation Blood Collection Tube



Retractable technology that withdraws needle into the tube using vacuum pressure




Activated Tube

The safety mechanism is an integral part of the device and **doesn't block the sight** of the needle tip.

Activation virtually eliminates exposure to both ends of the contaminated needle, **reducing the risk of needlestick injury**.

Allows users to **keep both hands behind the needle** when activating the needle retraction feature.

Instructions for Use ▶




Sol-Care™ Safety Blood Collection Tube Holder Set

Ref. Number	Box	Case
7001BCD	50	400

BLOOD COLLECTION



Sol-M™ Multi-sample Needle Sol-Care™ Safety Multi-sample Needle



Low-angle bevel design
Allows easy penetration.

Color-coded caps and hubs
Allow easy identification of needle gauge size.

Not made with natural rubber latex.

Sol-M™

Multi-sample Needle

Sol-M™ Multi-Sample Needle makes it possible to collect several samples with a single vein puncture.

- Marker on needle cap allows for bevel-up identification



Sol-M™ Multi-sample Needle

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Box Qty	Case Qty
MSN2310	23 G	0.6 mm	1"	25 mm	100	1000
MSN2210	22 G	0.7 mm	1"	25 mm	100	1000
MSN2215	22 G	0.7 mm	1 1/2"	38 mm	100	1000
MSN2110	21 G	0.8 mm	1"	25 mm	100	1000
MSN2115	21 G	0.8 mm	1 1/2"	38 mm	100	1000
MSN2010	20 G	0.9 mm	1"	25 mm	100	1000
MSN2015	20 G	0.9 mm	1 1/2"	38 mm	100	1000
VMSN1810	18 G	1.2 mm	1"	25 mm	100	1000
VMSN1815	18 G	1.2 mm	1 1/2"	38 mm	100	1000

18G needle has been designed for veterinary use; contains a component made of natural rubber latex.



Sol-Care™

Safety Multi-sample Needle

Safety-engineered tool for multi-sample venous blood collection. The device is designed with a safety mechanism covering the needle after use. In the activated position, the needle cover guards against accidental needlestick during normal handling and disposal.

It is also offered with a pre-attached holder.

Sol-Care™ Safety Multi-sample Needle

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Box Qty	Case Qty
SMSN22125	22 G	0.7 mm	1 1/4"	32 mm	48	480
SMSN2215	22 G	0.7 mm	1 1/2"	38 mm	48	480
SMSN21125	21 G	0.8 mm	1 1/4"	32 mm	48	480
SMSN2115	21 G	0.8 mm	1 1/2"	38 mm	48	480
SMSN20125	20 G	0.9 mm	1 1/4"	32 mm	48	480
SMSN2015	20 G	0.9 mm	1 1/2"	38 mm	48	480

Sol-Care™ Safety Multi-sample Needle with Pre-attached Holder


Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Box Qty	Case Qty
SMSNH22125	22 G	0.7 mm	1 1/4"	32 mm	48	480
SMSNH2215	22 G	0.7 mm	1 1/2"	38 mm	48	480
SMSNH21125	21 G	0.8 mm	1 1/4"	32 mm	48	480
SMSNH2115	21 G	0.8 mm	1 1/2"	38 mm	48	480
SMSNH20125	20 G	0.9 mm	1 1/4"	32 mm	48	480
SMSNH2015	20 G	0.9 mm	1 1/2"	38 mm	48	480



BLOOD COLLECTION



Sol-Care™ Safety Blood Collection Needle



Safety-engineered device intended for venous blood collection. The safety mechanism can be activated immediately after blood draw, helping to protect you against needlestick injury.

It is also offered with a pre-attached holder.

Sol-Care™

● Safety Blood Collection Needle



CE 0123

Sol-Care™ Safety Blood Collection Needle with Luer Adapter

Sol-Care™ Safety Blood Collection Needle with Pre-attached Holder



Sol-Care™ Safety Blood Collection Needle

with Luer Adapter

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Tubing Length (nches)	Tubing Length (mm)	Box Qty	Case Qty
110201030035	25 G	0.5 mm	3/4"	19 mm	12"	305 mm	50	200
110201030033	23 G	0.6 mm	3/4"	19 mm	12"	305 mm	50	200
110201030031	21 G	0.8 mm	3/4"	19 mm	12"	305 mm	50	200
110201030034	25 G	0.5 mm	3/4"	19 mm	7"	178 mm	50	200
110201030032	23 G	0.6 mm	3/4"	19 mm	7"	178 mm	50	200
110201030030	21 G	0.8 mm	3/4"	19 mm	7"	178 mm	50	200



Sol-Care™ Safety Blood Collection Needle


with pre-attached holder

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Tubing Length (nches)	Tubing Length (mm)	Box Qty	Case Qty
110201030010	25 G	0.5 mm	3/4"	19 mm	12"	305 mm	25	200
110201030006	23 G	0.6 mm	3/4"	19 mm	12"	305 mm	25	200
110201030002	21 G	0.8 mm	3/4"	19 mm	12"	305 mm	25	200
110201030012	25 G	0.5 mm	3/4"	19 mm	7"	178 mm	25	200
110201030008	23 G	0.6 mm	3/4"	19 mm	7"	178 mm	25	200
110201030004	21 G	0.8 mm	3/4"	19 mm	7"	178 mm	25	200

▶ **BLOOD COLLECTION**



Accessories



Group of accessories compatible with other products of the Sol-M blood collection system.

Sol-M™

● Luer Adapter



The Sol-M™ Luer Adapter is designed to work in a variety of functions by securely connecting to female Luer connections.

Ref. Number	Product description	Box Qty	Case Qty
LA21920	SOL-M™ Luer Adapter 20G	100	1000

Sol-M™

● Latex-Free Tourniquet



Pre cut, latex free elastic tourniquet strips are packaged in an easy dispensing box that is convenient to use and store.

Ref. Number	Product description	Box Qty	Case Qty
TNQ002	Sol-M™ Latex-Free Tourniquet 1" x 18"	25	2500

Sol-M™

● Holder



Sol-M™ Blood Collection Tube Holder

For blood specimen collection. Compatible with standard and safety blood collection needles.

Sol-M™ Blood Culture Holder

For blood culture collection and fits all standard evacuated blood culture bottles. Threads compatible with Sol-M Safety Blood Collection Needles.

SOL-M™ Tube Holder with Luer Lock Adapter

For collecting specimens from IV lines or for urine specimens. Standard female Luer fitting. Compatible with IV sets, ports, and tubing.

Sol-M™ Blood Transfer Holder

For the transfer of blood specimen from a syringe to vacuum blood tube. Female Luer Slip or Luer Lock compatible.

Ref. Number	Product description	Inner Bag Qty	Case Qty
110201050001	Sol-M™ Blood Collection Tube Holder	250	1000
110201050002	Sol-M™ Blood Culture Holder	40	320
110201050003	Sol-M™ Tube Holder with Luer Lock Adapter	50	800
BCTH002	Sol-M™ Blood Transfer Holder	50 in the box	800



Registered office

Sol-Millennium Europe Sp. z o.o.
Twarda 18, 00-105 Warsaw, Poland

Commercial office

Sol-Millennium Europe Sp. z o.o.
Warsaw HUB, Rondo Daszyńskiego 2B, 00-843 Warsaw, Poland

Registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIth Economic Department KRS 0000434783, NIP: 5252535720, EU VAT: PL5252535720, Share capital: 500.000 PLN

sol-m.com | info.eu@sol-m.com

EMEA-DIA-382-1022