



## Medecom MIP/MPR Mode

Available on Diagnostic workstation Med Mammo, Med Diag

## **Features**

- ✓ The MIP/MPR Mode is available on our diagnostic workstation Med Mammo & **Med Diag**
- ✓ MPR Mode is available from CT, MRI, ConeBeam exams
- ✓ Dedicated icon MPR gives instantly access to different functionalities
- ✓ Quick export of images of interest
- √ Very ergonomic mode

## **Functionalities**

- ✓ Easy access to double oblique mode to cut through volume following a different angle of incidence from standard views
- ✓ Point-to-point curve setting eases panoramic reconstruction, especially for dental applications
- ✓ Cut thickness and rake sizes (window of interest cuts) can simply be adjusted
- ✓ Smooth and intuitive:
  - **Navigation**
  - Reformatting
  - Exporting of interested images
  - Communication with PACS
- ✓ Axial, coronal, sagittal and 3D views simultaneously.

## **3D Navigation**

Different renders are available at your finger tips to visualize exams in a 3D volume form



MIP: Maximum intensity Projection MinIP: Minimum **Intensity Projection** AIP: Average **Intensity Projection** 



**Bones and Skin:** False color display brings out bones and skin



False color display brings out bones



False color display brings out the skin



Muscle: False color display brings out muscles

Medecom April 2020





**Medecom** is a private French company, created in 2000. Medecom counts more than 4 000 installations in 30 countries through a strong distribution network & OEM partnerships.